




To ensure that Innomed eNews is delivered to your inbox, please add eNews@innomed.net to your address book.
To view this page in your web browser on the internet, [click here](#).

APRIL 2017 – NEW PRODUCTS

Blair Acetabular Cup Positioner

Designed to help adjust the position of an acetabular cup

Designed by Christopher Blair, DO



Product No: 4159
Overall Length: 11.5" (29,2 cm)
Shaft Offset: 1" (2,54 cm)
Head Diameter: 1.75" Inside / 2" Outside (4,4 cm / 5,1cm)

Chung Weitlaner Retractor

Features long prongs allowing use of a small retractor in a small, but deep wound

Designed by Raymond Chung, MD

New 3x4 Prong Versions Available!



Particularly useful for volar approach to the distal radius in large patients and for mini-open rotator cuff surgery.

Product No's:



Prong lengths of 25mm and 30mm now available with either sharp or blunt tips, in both 3x2 or 3x4 prong versions.

BLUNT TIPS

5065 [2x3 Prongs]
Blade Depth: 25 mm
Overall Length: 4.5" (11,4 cm)

5065-01 [3x4 Prongs] **NEW VERSION!**
Blade Depth: 25 mm
Overall Length: 4.5" (11,4 cm)

5067 [2x3 Prongs]
Blade Depth: 30 mm
Overall Length: 4.5" (11,4 cm)

5067-01 [3x4 Prongs] **NEW VERSION**
Blade Depth: 30 mm
Overall Length: 4.5" (11,4 cm)

SHARPS TIPS

5066 [2x3 Prongs]
Blade Depth: 25 mm
Overall Length: 4.5" (11,4 cm)

5066-01 [3x4 Prongs] **NEW VERSION**
Blade Depth: 25 mm
Overall Length: 4.5" (11,4 cm)


5068 [2x3 Prongs]
Blade Depth: 30 mm
Overall Length: 4.5" (11,4 cm)


5068-01 [3x4 Prongs] **NEW VERSION**
Blade Depth: 30 mm
Overall Length: 4.5" (11,4 cm)

Doroodchi Coated Femoral Neck Mating Guide

Designed for controlled manipulation of femoral head/neck mating in SuperPATH THA approach

Designed by Hamidreza Doroodchi, MD





Product No: 3419
 Overall Length: 11.75 (29,8 cm)
 Blade Width: 1.125" (29 mm)

Goytia Osteotome Punch Tamp Assembly

Designed for removing a tibial bone plug to use as autograft for the femoral intramedullary alignment hole in total knee replacement

Designed by Robin Goytia, MD




Product No's:

5339-00 [Punch & Tamp Set]

Set Includes / Available Individually:

5339-01 [Osteotome Punch]
 Overall Length: 7.75" (19,7 cm)
 Outside Diameter: 16 mm
 Inside Diameter: 13.7 mm


5339-02 [Tamp]
 Overall Length: 7.75" (19,7 cm)
 Diameter: 12.3 mm



Kodkani Tissue Elevator Suture/Graft Passer

Designed for MPFL reconstruction basket weave technique, and helpful for mini-open ligament reconstruction surgeries for graft passage

Designed by Pranjal Kodkani, MD



Can also be used for:

Periosteum/soft tissue elevator or freer

Percutaneous passage of tendon/ligament graft/suture

Stripping tendon grafts off muscle

General orthopedics – periosteum elevator and spike

Product No: 1114
 Overall Length: 9.75" (24,8 cm)
 Handle Length: 4.25" (10,8 cm)
 Suture Hole: 2,5 mm x 13 mm

Light Source Cable Adapters

Designed to convert from our standard ACMI to either Storz or Wolf




Product No's:

8009-S ACMI to Storz Adapter
 8009-W ACMI to Wolf Adapter

Periarticular Reduction Forceps

Designed for reduction of intraarticular and periarticular fractures



Pointed ball tips help provide a secure hold in the bone despite minimal contact. Two sizes available.

Product No's:

1856 [Medium]

Overall Length: 14" (35,6 cm)

Width @ Tips Parallel: 10.5" (26,7 cm)

Maximum Jaw Opening @ Tips: 5.2" (13,2 cm)

1857 [Large]

Overall Length: 18.8" (47,8 cm)

Width @ Tips Parallel: 12" (30,5 cm)

Maximum Jaw Opening @ Tips: 7.3" (18,5 cm)

Taylor Retractors

**Product No's:****6330-01 [Standard]**

*Overall Length: 8" (20,3 cm)
Depth from Bend: 4" (10,2 cm)
Blade Width: 32 mm*

6330-02 [Deep]

*Overall Length: 9" (23 cm)
Depth from Bend: 5,5" (14 cm)
Blade Width: 32 mm*

6330-03 [Deep with Pin Guides]

*Overall Length: 9" (23 cm)
Depth from Bend: 5,5" (14 cm)
Blade Width: 32 mm
Guide for Pins Up To: .15" (3,5 mm)*

INSTRUMENTS FOR ORTHOPEDIC SURGERY

© 2017 Innomed, Inc., All rights reserved.